

Prime Priority Technology Ltd.

SpO2

Nellcor	Masimo	Philips	Mindray	Edan	Nihon Kohden
<p>TR-SP-TSR-0201</p> <p>Nellcor® Oximax: All Nellcor modulated patient</p>  <p>Nellcor Oximax reuse 9 pin adult finger clip sensor</p>	<p>TR-SP-TSR-0203</p> <p>Nellcor® Oximax: All Nellcor modulated patient</p>  <p>Nellcor Oximax reuse 9 pin adult finger silicone sensor</p>	<p>TR-SP-TSD-0201</p> <p>Nellcor® Oximax: All Nellcor modulated patient</p>  <p>Nellcor Oximax disposable 9 pin adult finger SpO2 sensor</p>	<p>TR-SP-TSD-0204</p> <p>Nellcor® Oximax: All Nellcor modulated patient</p>  <p>Nellcor Oximax disposable 9 pin neonate SpO2 sensor</p>		
<p>TR-SP-TSD-0901</p> <p>Masimo M-LNCS modulated</p>  <p>Masimo M-LNCS disposable SpO2 sensor</p>	<p>TR-SP-TSR-0901</p> <p>Masimo M-LNCS modulated</p>  <p>Masimo M-LNCS reuse SpO2 sensor</p>	<p>TR-SP-TSR-6403</p> <p>Masimo M-LNCS modulated</p>  <p>Masimo M-LNCS SpO2 connection cable</p>	<p>TR-SP-TSR-0401</p> <p>Masimo LNCS modulated</p>  <p>Masimo LNCS reuse adult finger SpO2 sensor</p>	<p>TR-SP-TSR-3002</p> <p>Masimo LNCS modulated</p>  <p>Masimo LNCS SpO2 connection cable</p>	
<p>TR-SP-TSR-4701-3</p> <p>Philips : M, MP, VM series Defibrillator Monitor M4735A</p>  <p>Philips SpO2 direct probe adult finger clip</p>	<p>TR-SP-TSR-4703-3</p> <p>Philips : M, MP, VM series Defibrillator Monitor M4735A</p>  <p>Philips SpO2 direct probe adult finger silicone</p>	<p>TR-SP-TSR-4705-3</p> <p>Philips : M, MP, VM series Defibrillator Monitor M4735A</p>  <p>Philips SpO2 direct probe adult ear clip</p>	<p>TR-SP-TSR-4702</p> <p>Philips : M, MP, VM series Defibrillator Monitor M4735A</p>  <p>Philips SpO2 D shape 8 pin connection cable</p>		
<p>TR-SP-TSR-3601-3</p> <p>Mindray: Beneview T5 iPM 8, iPM 10, iPM 12, VS900</p>  <p>Mindray SpO2 direct probe</p>	<p>TR-SP-TSR-3603</p> <p>Mindray: Beneview T5 iPM 8, iPM 10, iPM 12, VS900</p>  <p>Mindray SpO2 connection cable</p>	<p>TR-SP-TSR-3402X(80°)</p> <p>Mindray: PM8000, VS800 (Oximax module)</p>  <p>Mindray SpO2 connect cable 80° (Oximax module)</p>	<p>TR-SP-TSR-3402</p> <p>Mindray: PM800, VS800 (Mindray module)</p>  <p>Mindray SpO2 direct cable (Mindray module)</p>		
<p>TR-SP-TSR-6401-3</p> <p>Mindray: PM8000, VS800 Edan : M3</p>  <p>SpO2 direct probe (Nellcor)</p>	<p>TR-SP-TSR-6405-3</p> <p>Mindray: PM8000, VS800 Edan : M3</p>  <p>SpO2 direct cable ear clip</p>	<p>TR-SP-TSR-4201-3</p> <p>Nihon Kohden : BSM series</p>  <p>Nihon Kohden SpO2 direct probe, adult finger clip</p>	<p>TR-SP-TSR-4205-3</p> <p>Nihon Kohden : BSM series</p>  <p>Nihon Kohden SpO2 direct probe, adult ear clip</p>		

NIBP

TR-BP-TSR-1023S-01 (Standard)

Size (cm)

Neo : 10-15
Pedi : 14-21.5

Adult
Small : 20.5-28
Standard : 27-35
Large : 34-43
X-Large : 42-54



27-35cm adult reuse BP cuff, single hose

Single hose BP

TR-BP-TSD-1023S-01



27-35cm adult disposal BP cuff, single hose

Mindray, Edan, HF

TR-BP-TBR-1023S-08



BP connection tube

Philips, Mindray, COMEN

TR-BP-TBR-4303



BP connection tube

Double hose BP

TR-BP-TBR-4203




Dinamap BP connection tube

ECG cable

Philips

TR-EG-E2900CI15



Philips: M, MP, VM series
Defibrillator Monitor M4735A

Philips 5 lead ECG direct cable, IEC

Mindray

TR-EG-E2200CI13



Mindray: Beneview T5, iPM 8, iPM 10, iPM 12, VS900

Mindray 3 lead ECG direct cable, IEC


TR-EG-E11CI5-85



Philips: M, MP, VM series
Defibrillator Monitor M4735A

Philips 5 lead ECG cable (patient leadwire), IEC clip

TR-EG-E2911U15



Philips: M, MP, VM series
Defibrillator Monitor M4735A

Philips 5 lead ECG connection cable, IEC

TR-EG-E05CI3-85



Mindray: Beneview T5, iPM 8, iPM 10, iPM 12, VS900

Mindray 3 lead ECG cable (patient leadwire), IEC clip

TR-EG-E2205U13



Mindray: Beneview T5, iPM 8, iPM 10, iPM 12, VS900

Mindray 3 lead ECG connection cable, IEC

12 Lead

TR-EG-E0062BI03



Philips, Mindray, COMEN

Edan 12 lead ECG direct cable, IEC banana

TR-EG-E0069CI10-85



Philips, Mindray, COMEN

10 lead ECG cable (Half cable)

TR-EG-TS-01



Alligator clip for 4mm for banana pin

Philips Pogewriter Trim series
12 Lead Monitoring ECG Machine


TR-EG-E0064BI03



12 lead ECG direct cable for Philips TRIM series, IEC banana

TEMP

TR-TP-TP-01



All 2 pin temperature module

Disposable 2-pin rectal temp sensor

TR-TP-TP-02



All 2 pin temperature module

Disposable 2-pin skin temp sensor

TR-TP-TPR-47



Philips: M, MP, VM series
Defibrillator Monitor M4735A

Mindray/Philips 2-pin Temp adapter cable, Philips M1029A Mindray T5/T8 beneview

Mindray: Beneview T5, iPM 8, iPM 10, iPM 12, VS900

TR-TP-TPR-03



Mindray 2-pin Temp adapter cable, Mindray PM7000,8000,9000

Etco2

NT-CO-PH-CA01



Mindray: PM8000, VS800
Edan: M3

Respironics CARNOSTAT Mainstream CO2 Module

NT-CO-PH-CA02



Philips: M, MP, VM series
Defibrillator Monitor M4735A

Philips compatible external Mainstream CO2 sensor

NT-CO-PH-AD01



Respironics Adult Mainstream Airway Adaptor