

PADECG

PC ECG

Versatile

Mobile ECG Solution



EDAN

A world of potential

PADECG PC ECG

PADECG consists of a professional ECG software and a wireless ECG acquisition module. It provides a mobile solution to meet the demands of caregivers in a mobile treatment environment. EDAN's reliable SEMIP® algorithm is also provided in it, which allows you to have a fast and confident ECG assessment before sending the ECG report to SE-1515 Data Management System for further analysis.

PADECG Paperless Workflow



- User-friendly Interfaces
- Reliable Interpretation
- Extensive Communication

Touch Screen

- Operating PADECG is just as simple as using your smart phone. Convenient gestures such as long pressing, un-pinching, and dragging are all supported to enable users a better experience.

ECG Sampling

- Pre-sampling (only for Android) and real-time sampling modes are offered to cater to different situation, no matter it is an emergent case or general ward rounds.


Filter

- Armed with the comprehensive filters, you can acquire superb ECG traces with full confidence just within several touches.

Algorithm

- With EDAN's professional ECG algorithm SEMIP®, PADECG gives accurate diagnosis results and offers doctors a reliable reference.

★ *Tablet is not included in PADECG. Customers need to purchase tablet by themselves. The recommended tablets are: Samsung Galaxy Tab 3, Samsung Galaxy Tab PRO; iPad Series & iPad mini Series (iOS 7.0 and above).*



Prime Priority Technology Limited

Tel: **3703 9600**

Fax: 3011 5084

Email: prime.priority.tech@gmail.com



CE 0123

ENG-ECG-PADECG-V1.4-20180125