



Prime Priority Technology Limited

U 80 R

Large LCD size: 80 x 58mm with digit height 25mm
Product size: 139x95x34mm

Acrylic front panel and Ultra-thin cosmetic design

- 2x90 memory
- WHO indicator
- Year/month/date/time function
- Accuracy: Pressure +/-3mmHg; Pulse +/-5%
- Automatic power-off function
- Power by 4 AAA size alkaline batteries or AC/DC power supply
- 3 times average
- IHB indicator



Address: Room 1306, 13/F, Shing Yip Industrial Building, 19-21 Shing Yip Street, Kwun Tong, H.K.

Email: prime.priority.tech@gmail.com

Tel: 3703 9600 Fax: 3011 5084