



VeinNavi™ Simple Visibility

VeinNavi NAVI-60 vein finder provides a real time accurate image of the vasculature on the surface of the skin, which helps clinicians to verify the status of the vein and reduces the risk of venipuncture.

With the help of NAVI-60 vein finder, hospital may gain the benefit of: reduce costs for healthcare, improve infection control and improve patient safety and satisfaction.

Portable:

Hand hold type veinfinder, easy to carry, non-contact to patient.

Colorful:

Great image with 5 different colors optional, suitable for different types of skin.

Depth indicator:

The depth of the vein will be detected and shown on the side.



Specifications

Resolution	856*480 pixel
Image Light	10-25 Lumens 4 level adjustable
Color	5 colors
Image Size	3 Sizes
Vein depth detection	Yes
Screen	Yes
Charge during operation	Yes
Size	225*68*62mm
Weight	500g

+ APPLICATION



Plastic Surgery

Sclerotherapy for treating
spider telangiectasia
Cosmetic injections on the face



Vascular Surgery

Sclerotherapy for treating
spider telangiectasia
Vein stripping treating



Radiology

Limited viable patient access
Should puncture the big
and straight veins



Oncology

Fewer options for successful
access because of repeated
infusions and procedures
and medication



Pediatrics

Small veins
"Baby fat"
Anxiety of the children
and their families.

Prime Priority Technology Limited

Room 1306,13/F, Shing Yip Industrial Building, 19-21 Shing Yip Street, Kwun Tong, Hong Kong

Tel: 3703 9600

Fax: 3011 5084

Email: prime.priority.tech@gmail.com

