

yuwell

INFRARED EAR THERMOMETER

YHT107 



PRODUCT PARAMETERS

- Body measurement range : 34°C~42.2°C
- Temperature unit: °C/°F
- Time interval of each measurement: ≤12s
- Automatic power off time: 60s±10s
- Measurement duration: ≤2s
- Number of earmuffs: 20
- Memory retrieval: 10 memory groups

PACKAGING PARAMETERS

QTN: 40PCS/CTN
G.W./N.W: 10/3kg
Size: 67×37.5×26.5cm

FEVER GUIDANCE

DI-DI DI DI 



FAST MEASUREMENT

1s

BLUETOOTH (ONLY YHT107)



SPECIALIZED EARMUFFS



PRIME PRIORITY TECHNOLOGY LIMITED

Tel : (852) 3703 9600

Email : prime.priority.tech@gmail.com

www.primeprioritytech.com